

For Research Use Only

ZDHHC7 Polyclonal antibody

Catalog Number: 25985-1-AP



Basic Information

Catalog Number:

25985-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23185

GenBank Accession Number:

BC018772

GeneID (NCBI):

55625

UNIPROT ID:

Q9NXF8

Full Name:

zinc finger, DHHC-type containing 7

Calculated MW:

345 aa, 39 kDa

Observed MW:

45-49 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PC-3 cells, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-3 cells were subjected to SDS PAGE followed by western blot with 25985-1-AP (ZDHHC7 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.