

For Research Use Only

MTMR7 Polyclonal antibody

Catalog Number: 25990-1-AP



Basic Information

Catalog Number:

25990-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23170

GenBank Accession Number:

BC130577

GeneID (NCBI):

9108

UNIPROT ID:

Q9Y216

Full Name:

myotubularin related protein 7

Observed MW:

85 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:16000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Neuro-2a cells, mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25990-1-AP (MTMR7 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.