

For Research Use Only

# MTMR7 Polyclonal antibody

Catalog Number: 25990-1-AP



## Basic Information

<b>Catalog Number:</b> 25990-1-AP	<b>GenBank Accession Number:</b> BC130577	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 800 µg/ml	<b>GeneID (NCBI):</b> 9108	<b>Recommended Dilutions:</b> WB 1:2000-1:16000
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9Y216	
<b>Isotype:</b> IgG	<b>Full Name:</b> myotubularin related protein 7	
<b>Immunogen Catalog Number:</b> AG23170	<b>Observed MW:</b> 85 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : Neuro-2a cells, mouse brain tissue, rat brain tissue
<b>Species Specificity:</b> human, mouse, rat	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 25990-1-AP (MTMR7 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.