

For Research Use Only

GLT25D2 Polyclonal antibody

Catalog Number: 25993-1-AP **1 Publications**



Basic Information

Catalog Number:

25993-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23360

GenBank Accession Number:

BC035672

GeneID (NCBI):

23127

UNIPROT ID:

Q8IYK4

Full Name:

glycosyltransferase 25 domain
containing 2

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB: SH-SY5Y cells, Neuro-2a cells

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|------------------|-----------|----------------|-------------|
| Krista A Geister | 31101663 | Dis Model Mech | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 25993-1-AP (GLT25D2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.