

For Research Use Only

SLIRP Polyclonal antibody, PBS Only

Catalog Number: 26006-1-PBS

Featured Product



Basic Information

Catalog Number:

26006-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22945

GenBank Accession Number:

BC017895

GeneID (NCBI):

81892

UNIPROT ID:

Q9GZT3

Full Name:

chromosome 14 open reading frame

156

Calculated MW:

109 aa, 12 kDa

Observed MW:

12 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Background Information

SLIRP, also named as C14orf156, DC23, DC50 and PD04872. It is an RNA-binding protein on steroid receptor RNA activator (SRA), and its expression is ubiquitous in human normal tissues. SLIRP can affect mitochondrial mRNA transcription and energy conversion. In addition, the loss or weakening of SLIRP destroys part of the structure of mitochondria, affects mitochondria energy metabolism, and increases sperm or cell oxidative stress damage (PMID: 33150185). SLIRP has 2 isoforms with the molecular mass of 12 kDa.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 26006-1-AP (SLIRP antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26006-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 26006-1-AP (SLIRP antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26006-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using 26006-1-AP (SLIRP antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26006-1-PBS in a different storage buffer formulation.




Immunofluorescent analysis of (4% PFA) fixed PC-3 cells using SLIRP antibody (26006-1-AP) at dilution of 1:200 and CoralLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 26006-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed PC-3 cells using SLIRP antibody (26006-1-AP) at dilution of 1:200 and CoralLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 26006-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using 26006-1-AP (SLIRP antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26006-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 26006-1-AP (SLIRP antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26006-1-PBS in a different storage buffer formulation.