

For Research Use Only

FAM116B Polyclonal antibody

Catalog Number: 26024-1-AP



Basic Information

Catalog Number:

26024-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22053

GenBank Accession Number:

BC031069

GeneID (NCBI):

414918

UNIPROT ID:

Q8NEG7

Full Name:

family with sequence similarity 116,
member B

Calculated MW:

585 aa, 66 kDa

Observed MW:

66 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, A549 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 26024-1-AP (FAM116B Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.