

For Research Use Only

CEP19 Polyclonal antibody

Catalog Number: 26036-1-AP **5 Publications**



Basic Information

Catalog Number: 26036-1-AP	GenBank Accession Number: BC007827	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 84984	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q96LKO	
Isotype: IgG	Full Name: chromosome 3 open reading frame 34	
Immunogen Catalog Number: AG22035	Calculated MW: 163 aa, 19 kDa	
	Observed MW: 23 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse kidney tissue, mouse liver tissue
Cited Applications: WB, IP, IF	
Species Specificity: human, mouse, rat	
Cited Species: human, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Quynh P H Nguyen	33038334	Dev Cell	IF
Maria Chatzifrangkeskou	36347932	Sci Rep	IF
Tomoharu Kanie	28625565	Dev Cell	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 26036-1-AP (CEP19 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.