

For Research Use Only

CEP19 Polyclonal antibody

Catalog Number: 26036-1-AP

5 Publications



Basic Information

Catalog Number:

26036-1-AP

Concentration:

700 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22035

GenBank Accession Number:

BC007827

GeneID (NCBI):

84984

UNIPROT ID:

Q96LKO

Full Name:

chromosome 3 open reading frame 34

Calculated MW:

163 aa, 19 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : mouse kidney tissue, mouse liver tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Quynh P H Nguyen	33038334	Dev Cell	IF
Maria Chatzifrangkeskou	36347932	Sci Rep	IF
Tomoharu Kanie	28625565	Dev Cell	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 26036-1-AP (CEP19 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.