

For Research Use Only

GCKR Polyclonal antibody, PBS Only

Catalog Number: 26042-1-PBS



Basic Information

Catalog Number:

26042-1-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22033

GenBank Accession Number:

BC130481

GeneID (NCBI):

2646

UNIPROT ID:

Q14397

Full Name:

glucokinase (hexokinase 4) regulator

Calculated MW:

625 aa, 69 kDa

Observed MW:

68 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26042-1-AP (GCKR antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26042-1-PBS in a different storage buffer formulation.