

For Research Use Only

# IL-6 Polyclonal antibody

Catalog Number: 26049-1-AP

5 Publications



## Basic Information

Catalog Number:

26049-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23467

GenBank Accession Number:

NM\_031168

GeneID (NCBI):

16193

UNIPROT ID:

P08505

Full Name:

interleukin 6

Calculated MW:

24 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC, ELISA

Species Specificity:

human

Cited Species:

human, mouse, bovine

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Zhigang Liu	36288667	Bioorg Chem	WB
Zhao Chen	36410585	Free Radic Biol Med	ELISA
Jiaxuan Zou	40250625	J Control Release	IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

