

For Research Use Only

SLFN11 Polyclonal antibody, PBS Only

Catalog Number: 26060-1-PBS



Basic Information

Catalog Number:

26060-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22030

GenBank Accession Number:

BC052586

GeneID (NCBI):

91607

UNIPROT ID:

Q7Z7L1

Full Name:

schlafen family member 11

Calculated MW:

901 aa, 103 kDa

Observed MW:

103 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 26060-1-AP (SLFN11 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26060-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A431 cells using 26060-1-AP (SLFN11 antibody), at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 26060-1-PBS in a different storage buffer formulation.



IP result of anti-SLFN11 (IP:26060-1-AP, 4ug; Detection:26060-1-AP 1:800) with HEK-293 cells lysate 1610 ug. This data was developed using the same antibody clone with 26060-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human ovary cancer tissue slide using 26060-1-AP (SLFN11 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26060-1-PBS in a different storage buffer formulation.