

For Research Use Only

# UBR1 Polyclonal antibody, PBS Only

Catalog Number:26069-1-PBS



## Basic Information

Catalog Number:

26069-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22076

GenBank Accession Number:

BC113505

GeneID (NCBI):

197131

UNIPROT ID:

Q8IWW7

Full Name:

ubiquitin protein ligase E3  
component n-recognin 1

Calculated MW:

1749 aa, 200 kDa

Observed MW:

200 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26069-1-AP (UBR1 antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26069-1-PBS in a different storage buffer formulation.



IP result of anti-UBR1 (IP:26069-1-AP, 4ug; Detection:26069-1-AP 1:300) with mouse skeletal muscle tissue lysate 4000ug. This data was developed using the same antibody clone with 26069-1-PBS in a different storage buffer formulation.