

For Research Use Only

RLN3 Polyclonal antibody

Catalog Number: 26075-1-AP **8 Publications**



Basic Information

Catalog Number: 26075-1-AP	GenBank Accession Number: BC140935	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 117579	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q8WXF3	
Isotype: IgG	Full Name: relaxin 3	
Immunogen Catalog Number: AG23382		

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human stomach tissue,
Cited Applications: WB	
Species Specificity: human, mouse	
Cited Species: human, mouse	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Dongyuan Wang	36218134	Nanoscale	WB
Jinhong Zhu	36621836	Cancer Med	WB
Yulian Chen	35027826	Int J Nanomedicine	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 26075-1-AP (RLN3 antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 26075-1-AP (RLN3 antibody) at dilution of 1:200 (under 40x lens).