

For Research Use Only

RLN3 Polyclonal antibody

Catalog Number: 26075-1-AP

9 Publications



Basic Information

Catalog Number:

26075-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23382

GenBank Accession Number:

BC140935

GeneID (NCBI):

117579

UNIPROT ID:

Q8WXF3

Full Name:

relaxin 3

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

Applications

Tested Applications:

IHC, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human stomach tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Dongyuan Wang	36218134	Nanoscale	WB
Jinhong Zhu	36621836	Cancer Med	WB
Yulian Chen	35027826	Int J Nanomedicine	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 26075-1-AP (RLN3 antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human stomach tissue slide using 26075-1-AP (RLN3 antibody) at dilution of 1:200 (under 40x lens).