

For Research Use Only

C19orf20 Polyclonal antibody

Catalog Number: 26113-1-AP



Basic Information

Catalog Number:

26113-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23546

GenBank Accession Number:

BC009520

GeneID (NCBI):

91978

UNIPROT ID:

Q6ZTW0

Full Name:

chromosome 19 open reading frame

20

Observed MW:

29 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : human plasma tissue, mouse heart tissue, mouse thymus tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human plasma were subjected to SDS PAGE followed by western blot with 26113-1-AP (C19orf20 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.