

For Research Use Only

C1orf103 Polyclonal antibody

Catalog Number: 26115-1-AP

Featured Product

3 Publications



Basic Information

Catalog Number:

26115-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23494

GenBank Accession Number:

BC008115

GeneID (NCBI):

55791

UNIPROT ID:

Q5T3J3

Full Name:

chromosome 1 open reading frame

103

Observed MW:

85 kDa, 27 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : L02 cells, mouse heart tissue, mouse liver tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Chenlu Wang	38976714	PLoS Pathog	WB
Darina Šikrová	37380887	Commun Biol	WB
Kohei Hamanaka	32467133	Neurology	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



L02 cells were subjected to SDS PAGE followed by western blot with 26115-1-AP (C1orf103 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.