

For Research Use Only

DMRTB1 Polyclonal antibody

Catalog Number: 26118-1-AP



Basic Information

Catalog Number:

26118-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23639

GenBank Accession Number:

BC029566

GeneID (NCBI):

63948

UNIPROT ID:

Q96MA1

Full Name:

DMRT-like family B with proline-rich
C-terminal, 1

Observed MW:

40-45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human testis tissue,

Background Information

DMRTB1, also named as Doublesex- and mab-3-related transcription factor B1, is a 342 amino acid protein, which contains 1 DM DNA-binding domain and belongs to the DMRT family. DMRTB1 is expressed in testis.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 26118-1-AP (DMRTB1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.