

For Research Use Only

LSM11 Polyclonal antibody

Catalog Number: 26119-1-AP **1 Publications**



Basic Information

Catalog Number: 26119-1-AP	GenBank Accession Number: BC126449	Purification Method: Antigen affinity purification
Size: 350 µg/ml	GeneID (NCBI): 134353	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: P83369	
Isotype: IgG	Full Name: LSM11, U7 small nuclear RNA associated	
Immunogen Catalog Number: AG23281	Observed MW: 45-50 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2 cells
Cited Applications: WB	
Species Specificity: human	
Cited Species: human, mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Francesco S Ielasi	36447390	Nucleic Acids Res	WB


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various cells were subjected to SDS PAGE followed by western blot with 26119-1-AP (LSM11 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.