

For Research Use Only

LSM11 Polyclonal antibody

Catalog Number: 26119-1-AP

1 Publications



Basic Information

Catalog Number:

26119-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23281

GenBank Accession Number:

BC126449

GeneID (NCBI):

134353

UNIPROT ID:

P83369

Full Name:

LSM11, U7 small nuclear RNA associated

Observed MW:

45-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB : HeLa cells, HepG2 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Francesco S Ielasi	36447390	Nucleic Acids Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various cells were subjected to SDS PAGE followed by western blot with 26119-1-AP (LSM11 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.