

For Research Use Only

ZFP90 Polyclonal antibody

Catalog Number: 26120-1-AP



Basic Information

Catalog Number:

26120-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23381

GenBank Accession Number:

BC140886

GeneID (NCBI):

146198

UNIPROT ID:

Q8TF47

Full Name:

zinc finger protein 90 homolog
(mouse)

Calculated MW:

73 kDa

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat, mouse

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 26120-1-AP (ZFP90 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.