

For Research Use Only

ECHDC2 Polyclonal antibody

Catalog Number: 26126-1-AP



Basic Information

Catalog Number: 26126-1-AP	GenBank Accession Number: BC044574	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 55268	Recommended Dilutions: WB 1:1000-1:6000
Source: Rabbit	UNIPROT ID: Q86YB7	
Isotype: IgG	Full Name: enoyl Coenzyme A hydratase domain containing 2	
Immunogen Catalog Number: AG23510	Observed MW: 29-31 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : rat kidney tissue,
Species Specificity: human, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat kidney tissue were subjected to SDS PAGE followed by western blot with 26126-1-AP (ECHDC2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.