

For Research Use Only

ECHDC2 Polyclonal antibody, PBS Only

Catalog Number: 26126-1-PBS



Basic Information

Catalog Number:

26126-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23510

GenBank Accession Number:

BC044574

GeneID (NCBI):

55268

UNIPROT ID:

Q86YB7

Full Name:

enoyl Coenzyme A hydratase domain containing 2

Observed MW:

29-31 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Background Information

ECHDC2 (Enoyl-CoA Hydratase Domain Containing 2) is a protein-coding gene that encodes a mitochondrial enzyme. ECHDC2 is involved in fatty acid β -oxidation and branched-chain amino acid catabolism.

Storage

Storage:

Store at -80°C .

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat kidney tissue were subjected to SDS PAGE followed by western blot with 26126-1-AP (ECHDC2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26126-1-PBS in a different storage buffer formulation.