

For Research Use Only

C14orf119 Polyclonal antibody

Catalog Number: 26141-1-AP



Basic Information

Catalog Number:

26141-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23415

GenBank Accession Number:

BC009645

GeneID (NCBI):

55017

UNIPROT ID:

Q9NWQ9

Full Name:

chromosome 14 open reading frame

119

Observed MW:

13-16 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : A431 cells, A549 cells, HepG2 cells, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26141-1-AP (C14orf119 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.