

For Research Use Only

C2orf60 Polyclonal antibody

Catalog Number: 26145-1-AP



Basic Information

Catalog Number: 26145-1-AP	GenBank Accession Number: BC136837	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 129450	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: A2RUC4	
Isotype: IgG	Full Name: chromosome 2 open reading frame 60	
Immunogen Catalog Number: AG23452	Observed MW: 40-45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB: Jurkat cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 26145-1-AP (C2orf60 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.