For Research Use Only

Cytokeratin 71 Polyclonal antibody

Catalog Number: 26149-1-AP



Basic Information

Catalog Number:

GenBank Accession Number:

Purification Method: Antigen affinity purification

WB: 1:500-1:2000

26149-1-AP

GeneID (NCBI):

Recommended Dilutions:

Concentration: 1000 ug/ml

112802

Source: Rabbit

UNIPROT ID: Q3SY84

BC103917

Isotype:

Full Name:

IgG

keratin 71

Immunogen Catalog Number:

Observed MW: 57 kDa

Applications

Tested Applications:

WB, ELISA

Positive Controls:
WB: mouse skin tissue,

Species Specificity:

human, mouse

Background Information

KRT71 belongs to a family of type II keratins that are specifically expressed in the inner root sheath of hair follicles. The type II keratins are clustered in a region of chromosome 12q13.

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 26149-1-AP (KRT71 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.