

For Research Use Only

RPS6KB2 Polyclonal antibody

Catalog Number: 26194-1-AP



Basic Information

Catalog Number:

26194-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23800

GenBank Accession Number:

BC000094

GeneID (NCBI):

6199

UNIPROT ID:

Q9UBS0

Full Name:

ribosomal protein S6 kinase, 70kDa, polypeptide 2

Calculated MW:

53 kDa

Observed MW:

54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MCF-7 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 26194-1-AP (RPS6KB2 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.