

For Research Use Only

RPS6KB2 Polyclonal antibody

Catalog Number: 26194-1-AP



Basic Information

Catalog Number: 26194-1-AP	GenBank Accession Number: BC000094	Purification Method: Antigen affinity purification
Size: 750 µg/ml	GeneID (NCBI): 6199	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q9UBS0	
Isotype: IgG	Full Name: ribosomal protein S6 kinase, 70kDa, polypeptide 2	
Immunogen Catalog Number: AG23800	Calculated MW: 53 kDa	
	Observed MW: 54 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : MCF-7 cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 26194-1-AP (RPS6KB2 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.