

For Research Use Only

# NT5DC3 Polyclonal antibody

Catalog Number: 26242-1-AP



## Basic Information

Catalog Number:

26242-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23962

GenBank Accession Number:

BC052293

GeneID (NCBI):

51559

UNIPROT ID:

Q86UY8

Full Name:

5'-nucleotidase domain containing 3

Observed MW:

63 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

## Applications

Tested Applications:

IP, ELISA

Species Specificity:

human, mouse

Positive Controls:

IP : mouse testis tissue, HepG2 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IP result of anti-NT5DC3 (IP:26242-1-AP, 4ug;  
Detection:26242-1-AP 1:300) with mouse testis  
tissue lysate 4000ug.