

For Research Use Only

PDE7A Polyclonal antibody

Catalog Number: 26277-1-AP **1 Publications**



Basic Information

| | | |
|---|--|--|
| Catalog Number: 26277-1-AP | GenBank Accession Number: BC126360 | Purification Method: Antigen affinity purification |
| Size: 300 µg/ml | GeneID (NCBI): 5150 | Recommended Dilutions: WB 1:500-1:2000 |
| Source: Rabbit | UNIPROT ID: Q13946 | |
| Isotype: IgG | Full Name: phosphodiesterase 7A | |
| Immunogen Catalog Number: AG24334 | Calculated MW: 482 aa, 56 kDa | |
| | Observed MW: 48 kDa | |

Applications

| | |
|--|---|
| Tested Applications: WB, ELISA | Positive Controls: WB : Jurkat cells, K-562 cells |
| Cited Applications: WB | |
| Species Specificity: human | |
| Cited Species: rat | |

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-----------------------|-----------|---------------|-------------|
| Jose A Morales-Garcia | 31473904 | Mol Neurobiol | WB |


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 26277-1-AP (PDE7A Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.