

For Research Use Only

PDE7A Polyclonal antibody

Catalog Number: 26277-1-AP

1 Publications



Basic Information

Catalog Number:

26277-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24334

GenBank Accession Number:

BC126360

GeneID (NCBI):

5150

UNIPROT ID:

Q13946

Full Name:

phosphodiesterase 7A

Calculated MW:

482 aa, 56 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

rat

Positive Controls:

WB : Jurkat cells, K-562 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jose A Morales-Garcia	31473904	Mol Neurobiol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 26277-1-AP (PDE7A Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.