

For Research Use Only

CBX7 Polyclonal antibody

Catalog Number: 26278-1-AP

Featured Product

10 Publications



Basic Information

Catalog Number:

26278-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24588

GenBank Accession Number:

BC051773

GeneID (NCBI):

23492

UNIPROT ID:

O95931

Full Name:

chromobox homolog 7

Calculated MW:

28 kDa

Observed MW:

22-36 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IP, IF, IHC

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : mouse liver tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jiannan Ren	36396821	Br J Cancer	WB,IHC,IP
Ge Gao	35693076	Am J Cancer Res	WB
Wentao Liu	35342353	Int J Biol Sci	WB,IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 26278-1-AP (CBX7 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.