

For Research Use Only

# CCDC135 Polyclonal antibody, PBS Only

Catalog Number: 26307-1-PBS



## Basic Information

Catalog Number:

26307-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23635

GenBank Accession Number:

BC036667

GeneID (NCBI):

84229

UNIPROT ID:

Q8IY82

Full Name:

coiled-coil domain containing 135

Observed MW:

103 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26307-1-AP (CCDC135 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26307-1-PBS in a different storage buffer formulation.