

For Research Use Only

C1orf71 Polyclonal antibody

Catalog Number: 26322-1-AP

1 Publications



Basic Information

Catalog Number:

26322-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23732

GenBank Accession Number:

BC010228

GeneID (NCBI):

163882

UNIPROT ID:

Q6PJW8

Full Name:

chromosome 1 open reading frame 71

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse cerebellum tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Haoyu Liu	35662693	Front Pharmacol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with 26322-1-AP (C1orf71 Antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.