

For Research Use Only

GPAT2 Polyclonal antibody

Catalog Number: 26328-1-AP



Basic Information

Catalog Number:

26328-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23677

GenBank Accession Number:

BC068596

GeneID (NCBI):

150763

UNIPROT ID:

Q6NUI2

Full Name:

hypothetical protein LOC150763

Observed MW:

88 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:2000-1:10000

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

WB, IP, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, rat testis tissue

IP : rat testis tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26328-1-AP (GPAT2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



IP result of anti-GPAT2 (IP:26328-1-AP, 4ug; Detection:26328-1-AP 1:1000) with rat testis tissue lysate 1800 ug.