

For Research Use Only

GTSF1 Polyclonal antibody

Catalog Number: 26340-1-AP



Basic Information

Catalog Number:

26340-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23647

GenBank Accession Number:

BC021179

GeneID (NCBI):

121355

UNIPROT ID:

Q8WW33

Full Name:

gametocyte specific factor 1

Observed MW:

19-23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : K-562 cells, mouse testis tissue, rat testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26340-1-AP (GTSF1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.