

For Research Use Only

# H1FOO Polyclonal antibody

Catalog Number: 26347-1-AP



## Basic Information

Catalog Number:

26347-1-AP

Size:

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23688

GenBank Accession Number:

BC047943

GeneID (NCBI):

132243

UNIPROT ID:

Q8I2A3

Full Name:

H1 histone family, member O, oocyte-specific

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

## Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP: HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

