

For Research Use Only

LRRC42 Polyclonal antibody, PBS Only

Catalog Number: 26348-1-PBS



Basic Information

Catalog Number:

26348-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23710

GenBank Accession Number:

BC013101

GeneID (NCBI):

115353

UNIPROT ID:

Q9Y546

Full Name:

leucine rich repeat containing 42

Observed MW:

45-48 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

LRRC 42 (Leucine-Rich Repeating-Containing Protein 42) is a human protein, which is encoded by LRRC 42 gene and contains leucine-rich repetitive sequences. This protein exists in human body and is considered as a potential biomarker in many tumor types, which is involved in the key cellular process of cancer development and may be related to immune cell infiltration.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26348-1-AP (LRRC42 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26348-1-PBS in a different storage buffer formulation.