

For Research Use Only

# FOPNL Polyclonal antibody

Catalog Number: 26409-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

26409-1-AP

Size:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23991

GenBank Accession Number:

BC022321

GeneID (NCBI):

123811

UNIPROT ID:

Q96NB1

Full Name:

chromosome 16 open reading frame

63

Observed MW:

22 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

IF, IHC

Species Specificity:

human

Cited Species:

zebrafish

Positive Controls:

WB : HEK-293 cells,

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Shanshan Xie	30475641	FASEB J	IF, IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 26409-1-AP (FOPNL Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.