

For Research Use Only

KCNH1 Polyclonal antibody

Catalog Number: 26426-1-AP



Basic Information

Catalog Number: 26426-1-AP	GenBank Accession Number: BC113709	Purification Method: Antigen affinity purification
Concentration: 550 ug/ml	GeneID (NCBI): 3756	Recommended Dilutions: WB: 1:5000-1:50000 IHC: 1:400-1:1600
Source: Rabbit	UNIPROT ID: O95259	
Isotype: IgG	Full Name: potassium voltage-gated channel, subfamily H (eag-related), member 1	
Immunogen Catalog Number: AG16647	Calculated MW: 989 aa, 111 kDa Observed MW: 111 kDa	

Applications

Tested Applications: WB, IHC, ELISA	Positive Controls:
Species Specificity: human, mouse, rat	WB : mouse brain tissue, rat brain tissue IHC : rat brain tissue, mouse brain tissue
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Potassium voltage-gated channel subfamily H member 1 (KCNH1, also known as Kv10.1 and EAG) a member of the EAG (ether-à-go-go) family, which is involved in intracellular signaling, cell proliferation, and tumorigenesis (PMID: 35639255). KCNH1 is mainly expressed in the adult central nervous system. It also plays a role in controlling K⁺ flux that regulates resting membrane potential in excitable and non-excitable cells and is activated at the onset of myoblast differentiation (PMID: 29333676; 27029528; 9738473).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 26426-1-AP (KCNH1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.




IP result of anti-KCNH1 (IP:26426-1-AP, 4ug; Detection:26426-1-AP 1:5000) with rat brain tissue lysate 1720 ug.




IP result of anti-KCNH1 (IP:26426-1-AP, 4ug; Detection:26426-1-AP 1:5000) with mouse brain tissue lysate 2120 ug.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 26426-1-AP (KCNH1 antibody) at dilution of 1:800 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded rat brain tissue slide using 26426-1-AP (KCNH1 antibody) at dilution of 1:800 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded rat brain tissue slide using 26426-1-AP (KCNH1 antibody) at dilution of 1:800 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).