

For Research Use Only

RASSF3 Polyclonal antibody, PBS Only

Catalog Number: 26428-1-PBS



Basic Information

Catalog Number:

26428-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15801

GenBank Accession Number:

BC100952

GeneID (NCBI):

283349

UNIPROT ID:

Q86WH2

Full Name:

Ras association (RalGDS/AF-6)
domain family member 3

Calculated MW:

238 aa, 28 kDa

Observed MW:

28-30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IF/ICC, ELISA

Species Specificity:

human

Background Information

RASSF3 (Ras Association Domain Family Member 3) is an emerging tumor suppressor protein and a key member of the RASSF family, characterized by a conserved Ras Association (RA) domain. RASSF3 acts as a critical node in multiple tumor-suppressive pathways. It directly stabilizes and activates the p53 protein, thereby inducing cell cycle arrest and apoptosis. RASSF3 also contributes to the maintenance of genomic integrity by participating in the DNA damage response. Furthermore, it modulates major signaling cascades such as the Hippo and NF- κ B pathways, influencing cell proliferation, survival, and inflammation.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



LNCaP cells were subjected to SDS PAGE followed by western blot with 26428-1-AP (RASSF3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26428-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using RASSF3 antibody (26428-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 26428-1-PBS in a different storage buffer formulation.