

For Research Use Only

FAM98A Polyclonal antibody, PBS Only

Catalog Number: 26432-1-PBS



Basic Information

Catalog Number:

26432-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23220

GenBank Accession Number:

BC136770

GeneID (NCBI):

25940

UNIPROT ID:

Q8NCA5

Full Name:

family with sequence similarity 98, member A

Observed MW:

50-55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26432-1-AP (FAM98A antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26432-1-PBS in a different storage buffer formulation.