

For Research Use Only

ZDHHC12 Polyclonal antibody

Catalog Number: 26433-1-AP



Basic Information

Catalog Number: 26433-1-AP	GenBank Accession Number: BC048251	Purification Method: Antigen affinity purification
Size: 1000 µg/ml	GeneID (NCBI): 84885	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: Q96GR4	
Isotype: IgG	Full Name: zinc finger, DHHC-type containing 12	
Immunogen Catalog Number: AG23634	Observed MW: 29 kDa, 32 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC: human kidney tissue, human thyroid cancer tissue

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 26433-1-AP (ZDHHC12 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).