

For Research Use Only

ZDHHC12 Polyclonal antibody

Catalog Number: 26433-1-AP

1 Publications



Basic Information

Catalog Number:

26433-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23634

GenBank Accession Number:

BC048251

GeneID (NCBI):

84885

UNIPROT ID:

Q96GR4

Full Name:

zinc finger, DHHC-type containing 12

Observed MW:

29 kDa, 32 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IHC 1:50-1:500

Applications

Tested Applications:

IHC, ELISA

Cited Applications:

IP

Species Specificity:

human, mouse

Cited Species:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC : human kidney tissue, human thyroid cancer tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shuai Yang	39225180	Autophagy	IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 26433-1-AP (ZDHHC12 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).