

For Research Use Only

# ZBTB2 Polyclonal antibody

Catalog Number: 26436-1-AP



## Basic Information

Catalog Number:

26436-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23897

GenBank Accession Number:

BC020172

GeneID (NCBI):

57621

UNIPROT ID:

Q8N680

Full Name:

zinc finger and BTB domain  
containing 2

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 µg for 1.0-3.0 mg of total  
protein lysate

## Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP: HEK-293T cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IP result of anti-ZBTB2 (IP:26436-1-AP, 4ug;  
Detection:26436-1-AP 1:500) with HEK-293T cells  
lysate 2040 ug.