

For Research Use Only

Nesprin 4 Polyclonal antibody

Catalog Number: 26467-1-AP



Basic Information

Catalog Number:

26467-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24107

GenBank Accession Number:

BC052573

GeneID (NCBI):

163183

UNIPROT ID:

Q8N205

Full Name:

chromosome 19 open reading frame

46

Observed MW:

44 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, L02 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 26467-1-AP (Nesprin 4) antibody at dilution of 1:300 incubated at room temperature for 1.5 hours.