

For Research Use Only

# FAM161B Polyclonal antibody

Catalog Number: 26492-1-AP



## Basic Information

Catalog Number:

26492-1-AP

Size:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24566

GenBank Accession Number:

BC053909

GeneID (NCBI):

145483

UNIPROT ID:

Q96MY7

Full Name:

family with sequence similarity 161,  
member B

Observed MW:

72-74 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : THP-1 cells, U-251 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26492-1-AP (FAM161B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.