

For Research Use Only

FAM161B Polyclonal antibody, PBS Only

Catalog Number: 26492-1-PBS



Basic Information

Catalog Number:

26492-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24566

GenBank Accession Number:

BC053909

GeneID (NCBI):

145483

UNIPROT ID:

Q96MY7

Full Name:

family with sequence similarity 161, member B

Observed MW:

72-74 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26492-1-AP (FAM161B antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26492-1-PBS in a different storage buffer formulation.