

For Research Use Only

ZDHHC1 Polyclonal antibody

Catalog Number: 26545-1-AP **2 Publications**



Basic Information

Catalog Number: 26545-1-AP	GenBank Accession Number: BC021908	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 29800	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q8WTX9	
Isotype: IgG	Full Name: zinc finger, DHHC-type containing 1	
Immunogen Catalog Number: AG24028	Observed MW: 48-55 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB: A549 cells,
Cited Applications: WB, IHC, CoIP	
Species Specificity: human	
Cited Species: human	

Background Information

Catalog#26545-1-AP recognises 48-55 kDa bands.

Notable Publications

Author	Pubmed ID	Journal	Application
Cheng-Gong Luo	36589680	Oxid Med Cell Longev	WB
Qun Zhang	39069526	Cancer Gene Ther	WB, IHC, CoIP

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 26545-1-AP (ZDHC1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.