

For Research Use Only

ELOVL7 Polyclonal antibody

Catalog Number: 26615-1-AP



Basic Information

Catalog Number:

26615-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25008

GenBank Accession Number:

BC130310

GeneID (NCBI):

79993

UNIPROT ID:

A1L3X0

Full Name:

ELOVL family member 7, elongation of long chain fatty acids (yeast)

Calculated MW:

281 aa, 33 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26615-1-AP (ELOVL7 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.