

For Research Use Only

TMEM138 Polyclonal antibody

Catalog Number: 26642-1-AP



Basic Information

Catalog Number:

26642-1-AP

Concentration:

1100 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24487

GenBank Accession Number:

BC005201

GeneID (NCBI):

51524

UNIPROT ID:

Q9NPIO

Full Name:

transmembrane protein 138

Calculated MW:

162 aa, 19 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:200-1:800

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : hTERT-RPE1 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Methanol) fixed hTERT-RPE1 cells using TMEM138 antibody (26642-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CoraLite®594 ARL13B antibody (CL594-17711, red).